



## PREPARING THE CHAIR FOR USE




**Safety Note:** Do not release safety catches or adjust handle with passenger seated.



# DESCENDING THE STAIRS

1



- A. Release brakes.
- B. Move to top of stairs on 4 wheels
- C. Tilt chair backwards and position square to stairs.

Tip: Use the balance point for ease at this stage.

Safety Note: Use additional operators for small spaces, uneven surfaces etc.


2



Proceed forward over top step and down with tracks resting on stair nosings.

Tip: Communicate with passenger at every stage for reassurance.

3



- A. Operator maintains upright posture.
- B. Apply pressure towards stairs to control descent.





C. Stand chair on all wheels to turn and for travel on level surfaces.

Tip: You can also turn the chair inclined at the balance point

Tip: To ease forward tilt, stop with the rear wheels approx. 8cm (3") above landing.

Safety Note: Apply brakes if leaving chair and passenger unattended on level surfaces.

# CARRY FACILITY

<p><b>A</b></p> <ol style="list-style-type: none"><li>1. Raise rear carry / push handle</li><li>2. Prepare for use as Versa</li></ol>  <p>Tip: Refer to instruction label on rear handle.</p>	<ol style="list-style-type: none"><li>4. With passenger seated fasten and adjust 4 point harness straps if required</li></ol> 
<p><b>C</b></p> <ol style="list-style-type: none"><li>5. Tilt chair backwards at bottom step to raise front handles</li><li>6. Using 2 or more operators as necessary, lift the chair clear of the stairs and climb</li></ol>  <p>Tip: Also for spiral stairs, over obstacles, rubble etc.</p>	<p>also...</p>  <p>Superb mobility on level surfaces with extended rear push handle.</p> <p>Tip: Rear push handle can be folded down for stair descent if required.</p>