



Health, Safety and Wellbeing-Laser Safety Guidelines

General guidelines for all users

Contents

LASER SAFETY PROTOCOL	. 5
LEARN IT	.7
Introduction	. 9
What the protocol does not cover	. 9
Lasers: the basics	10
Hazards	10
Eye hazards	11
Classification	15
Roles	18
Training and requirements	20
Vision health check	20
University Laser Competency Matrix (v0.1)	22
Recommended Minimum Training:	24
Records	25
Laser Laboratory Authorised Personnel Training and Access Register	25
Laser Laboratory Authorised Personnel Training and Access Register	
	27
PURCHASE IT	27 29
PURCHASE IT	27 29 29
PURCHASE IT	27 29 29 30
PURCHASE IT Z Safe equipment. Z Inventory Z Appendix 1 Template to send to the LSO with laser details Z	29 29 30 31
PURCHASE IT Safe equipment. Safe equipment. Inventory Inventory Safe equipment is a send to the LSO with laser details Purchase request. Safe equipment is a send to the LSO with laser details	 27 29 29 30 31 33
PURCHASE IT Z Safe equipment. Z Inventory Z Appendix 1 Template to send to the LSO with laser details Z Purchase request. Z USE IT Z	 27 29 29 30 31 33 35
PURCHASE IT Z Safe equipment. Z Inventory Z Appendix 1 Template to send to the LSO with laser details Z Purchase request. Z USE IT Z General rules Z	 27 29 30 31 33 35 37
PURCHASE IT Safe equipment. Safe equipment. Inventory Inventory Appendix 1 Template to send to the LSO with laser details Purchase request. Purchase request. USE IT General rules Personal Protective Equipment (PPE). Safe the second s	 27 29 29 30 31 33 35 37 38
PURCHASE IT Safe equipment. Safe equipment. Inventory Appendix 1 Template to send to the LSO with laser details Purchase request. Purchase request. Purchase request. General rules Personal Protective Equipment (PPE). Storage spaces at the point of entry Storage spaces at the point of entry	 27 29 30 31 33 35 37 38 39
PURCHASE IT Safe equipment. Safe equipment. Inventory Inventory Appendix 1 Template to send to the LSO with laser details Purchase request. Purchase request. USE IT General rules Personal Protective Equipment (PPE) Storage spaces at the point of entry Good example of laser eyewear protection Good example of laser eyewear protection	 27 29 30 31 33 35 37 38 39 40

2/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Introduction
Risk assessment
During operation
Records
Laser survey 42
BUILD IT/MAINTAIN IT43
How to set up a laser lab
Risk assessment
Example risk assessment
Repairs and maintenance 59
Records
Temporary laser controlled areas
Controls during servicing
SECURE IT61
Physical security and access63
Information security63
Why do we need information security?63
What information needs to be kept securely?
Where is the information held?63
Using shared databases
Review of information security64
MOVE IT/DISPOSE OF IT65
Transfer
Disposal67
MEASURE IT69
Internal verification
Annual laser laboratory inspection71
Annual verification of the restricted lasers register
Documentation in each laser laboratory

3/79


Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

University laser laboratory evaluation checklist	72
EMERGENCY	75
Emergency procedures	77
Incident reports	77
Appendix 1. How to find the Eye Clinic	

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

LASER SAFETY PROTOCOL



Local operational guidance

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System



Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

7/79

HSW009.1

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Introduction

Lasers are special items of equipment that pose unique hazards to their users and others who may be nearby, so we need high quality, easy-to-follow instructions on how to use them safely. The aim of this guidance is to help staff and students operating lasers to follow good practices working in safety and in accordance with AS NZS IEC 60825.1 and AS NZS IEC 60825.14.



These guidelines are part of the Laser Safety protocol as shown below:

Local operational guidance

What the protocol does not cover

Medical or cosmetic lasers

Note that the laser safety protocol **does not** cover medical or cosmetic lasers. The use of such lasers in a medical setting is extremely complex and users must strictly comply with all applicable regulations and guidance – seek further advice.

Laser displays and shows

Lasers used in displays and shows are regulated under a specific standard (AS/NZS-IEC-608025.3 – Safety of laser products – Part 3: Guidance for laser displays and shows) which provides guidance on the planning and design, set-up and conduct of laser displays and shows that make use of high power lasers.

The laser power needed to produce effective theatrical or artistic displays in large spaces such as theatres, arenas, or architectural sites is great enough to pose a severe accidental exposure hazard, even when personal exposure is very brief. For this reason, sub-clause

9/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021 4.1.5 of IEC/TR 60825-14 specifies that only laser products that are Class 1, Class 2 or visible-beam Class 3R should be used for demonstration, display or entertainment purposes in unsupervised areas. Laser products of other classes should only be permitted under carefully controlled conditions and under the control of a trained experienced operator.

The guidance provides recommendations for safety for those laser displays or demonstrations that are shows, artistic displays, advertising or light sculptures, or museum pieces used to demonstrate optical principles, etc. Laser products available for use in a domestic environment or for use by people who cannot be expected to have received a suitable level of training should be Class 1, Class 2 or visible beam Class 3R. Therefore, such equipment is outside the scope of this guidance.

Lasers: the basics

The word "laser" is an acronym of Light Amplification by the Stimulated Emission of Radiation. A laser produces an intense beam of light with the unique properties of coherency, collimation and monochromaticity.

It is very common for lasers to be used in modern manufactured household products such as Blu-ray players, DVD players and optical computer drives. In most cases they are safe for the eyes because they are enclosed inside a protective housing, but some exposed lasers (especially lasers in research laboratories and high-powered laser pointers) can be dangerous if they are not properly used or set up.

A laser consists in a *gain medium*, capable of amplifying the light passing through, and an external system that provides energy, called *pumping*.

Lasers can be defined according to the material of the gain medium and also by the power output.

The laser can operate in:

- Continuous wave (CW) mode if the power output is continuous over time
- Pulsed mode if the power is turned on and off

According to the application that the laser will be used for, you might use one or the other. Usually the power used by lasers varies from a few milliwatts to several hundred watts.

Hazards

The hazards related to the use of lasers are:

• *Direct hazards,* caused directly by exposure of the eye or skin to the laser beam.


10/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

- *Indirect hazards,* caused by the interactions of the laser beam with reflecting objects in the environment.
- *Non-beam hazards,* including electrical hazards from the equipment, fire, hazardous fumes if directing lasers on to certain materials, human factor (fatigue, risk perception), computer software and more.

Eye hazards

Looking directly into a laser beam is known as intra-beam viewing, and all of the laser energy will enter the eye. If the laser is shone into a mirror, or onto another type of *specular reflector* (shiny surface) such as a wrist watch, a piece of shiny plastic or a metal object, the reflected laser energy may still be hazardous if the beam enters the eye. Depending on its power, even looking at the *diffuse* laser dot projected onto an object may be harmful. In all cases, using an optical device such as a pair of binoculars, a telescope, or a single lens reflex camera will concentrate laser energy within the eye and worsen any potential injury.




Different Ways to View Laser Beams

The energy (light transmission) entering the eye depends on the wavelength. This also defines different types of potential injuries.


11/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021



Light transmission through the eye according to the wavelength

Lasers operating in wavelengths between 400nm to 1400nm are known to easily cause retinal damage and potentially blind a person. This is why laser exposure standards are most strict for these wavelengths.



Percentage of the energy (transmission, red curve) which can reach the retina at the back of the eye.

12/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

However, UV (180-400nm) and infrared B and C (1400nm-1mm) lasers of Class 3B and 4 are also restricted because they can damage the cornea producing photokeratitis, photochemical cataract and corneal burns.

The level of laser radiation to which a person might be exposed without hazardous effects or adverse biological changes to the eyes or skin is called the *Maximum permissible exposure (MPE)*. The *MPE* can be used when determining the risk associated with the laser.

Most light sources tend to scatter optical radiation in all directions, and very little energy will actually enter the eye. In addition, because of the way the eye processes light, that energy will be widely spread across the light sensing cells in the back of the eye (*the macula*) as an image. With a laser, all of the light is aligned into a tight stream of energy known as a collimated laser beam. If this beam enters the eye, it will be concentrated even further onto the foveal region of the macula as a very intense dot.

If the intensity of the dot exceeds the MPE, the delicate light-sensing cells may be burned and permanently damaged. Extremely powerful lasers may even produce small steam explosions within the eye and cause catastrophic internal eye trauma.


Invisible lasers are extra-hazardous as the light energy being projected into the eye can still inflict damage. Being invisible, the laser beam will not trigger a person's natural blink reflex as would a bright, visible laser.

As an indication of the energy levels we are talking about:

- The energy from a 1 milliwatt laser shone into the eye is 167 times that of looking directly at the sun.
- The energy from a 5 watt laser shone into the eye is 835,000 times that of looking directly at the sun.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System



Intensity of a 1mW laser compared to a 100W light bulb

Fact: Some of the Class 4 lasers used in University research are 5,500,000 (five and a half million) times more powerful than the safe maximum permissible exposure (MPE)limit, and can cause permanent damage to unprotected eyes that are within 23 kilometres! Some of our very powerful lasers can also cause serious skin burns.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Classification

Lasers are now classified according to the following international system:

Laser class	Remarks
Class 1	No precautions are required. The laser is enclosed within a housing and cannot be directly viewed, or it is of such low power that it can be directly viewed without the need for eye protection.
Class 1C	Cosmetic/medical lasers. Although the emitted laser radiation may be at Class 3R, 3B or 4 levels, ocular exposures are prevented by one or more engineering means. Eye damage may occur if the protective measures are defeated.
Class 1M	The laser is eye safe with normal vision (you can directly view a Class 1M laser without having to wear eye protection). Eye damage may occur if you look at the laser through magnifying optical equipment (binoculars, telescopes, etc.).
Class 2	The laser must operate only in visible wavelengths, and is eye safe only because of the natural aversion/blink response (you instinctively blink when a bright light is shone in the eyes). The laser is safe to be viewed through magnifying optical equipment.
Class 2M	The laser must operate only in visible wavelengths, and is eye safe only because of the natural aversion/blink response (you instinctively blink when a bright light is shone in the eyes). Eye damage may occur if you look at the laser through magnifying optical equipment (binoculars, telescopes, etc.).
Class 3R	Lasers operating in visible wavelengths are generally eye safe because of the natural aversion/blink response. Eye damage may occur if you look at the laser through magnifying optical equipment (binoculars, telescopes, etc.). Invisible wavelengths <i>are not</i> eye safe as the blinking response will not be triggered.
Class 3B	Eye damage is very likely to occur. Generally safe for skin.
Class 4	Eye damage is very likely to occur. Not safe for skin.

Green: unrestricted lasers

Red: restricted lasers

15/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

For the purposes of the University Laser Safety Standard, lasers are grouped as follows:

- Unrestricted lasers are lasers which, if used under normal operational conditions and without modification (e.g. using lenses to focus the beam) are considered to pose a low or negligible risk. For the purpose of this protocol, all lasers that are classified as Class 1, 1C, 1M, 2, 2M or 3R in accordance with AS/NZS IEC 60825.1 are deemed to be unrestricted lasers.
- *Restricted lasers* are lasers with such potential to cause harm that they need to be carefully controlled. For the purpose of this protocol, all high powered laser pointers and lasers that are classified as Class 3B or Class 4 in accordance with AS/NZS IEC 60825.1 are deemed to be restricted lasers.

All lasers must be clearly labelled with a sticker that shows its classification. **If a sticker is not present**, the Laser Safety Officer should be contacted to determine whether the laser is safe to use.


Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System





Is my laser restricted?





Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Roles

A person's role is an indication of their competency when interacting with restricted lasers. People operating unrestricted lasers (Classes 1, 1C, 1M, 2, 2M and 3R) must comply with the lasers operating instructions and must not use the laser to blind, dazzle, or distract other people.

Roles within a workplace where restricted lasers are used are: Laser User, Laser Operator, Laser Supervisor, Laser Laboratory Coordinator and Laser Safety Officer.

Laser User: A laser user is a person who has had only basic laser safety awareness training, and is not yet competent enough to be designated a laser operator. They may only use restricted lasers when a laser supervisor is present. Laser users are normally students or inexperienced staff.

Laser Operator: Laser operators are people who have demonstrated competence to a level where they can work with minimum or no supervision on the specific lasers they have been trained to use. Laser operators are normally experienced staff, such as technicians and some postgraduate students. Laser operators may be authorised on an "as required" basis to act as a responsible person/monitor to "buddy" other operators using restricted risk lasers.

Laser Supervisor: Laser supervisors are people who are deemed competent to instruct others on how to use lasers that they are authorised to use. They can supervise users (such as staff or students) who have not yet achieved the levels of competency required to be laser operators. They may also be required to administer day-to-day laser activities within a workplace.

Laser Laboratory Coordinator: Laser laboratory coordinators (LLC) are experienced laser supervisors who have the authority to deem people competent in laser use and operation. They are able to instruct others on how to use lasers that they are authorised to use. They can supervise users (such as staff or students) who have not yet achieved the levels of competency required to be laser operators.

Laser safety officer: The laser safety officer (LSO) is a competent person who is appointed by the Associate Director, Health, Safety and Wellbeing, on behalf of the University. This person is required to be knowledgeable in the evaluation and control of laser hazards, and is responsible for overseeing the control of laser hazards at the University. Duties of the LSO include, but are not limited to, the following:

 Approve laser laboratories and other workplaces where restricted lasers are being used

- Maintaining a register (inventory) of all restricted lasers at the University
- Approving the acquisition of any restricted laser for the University (regardless of whether the laser is purchased, hired, borrowed, gifted, manufactured, converted, etc.).
- Providing information and specialist advice on laser management and operation as requested or required.
- Conducting periodic inspections of laser laboratories and other workplaces where restricted lasers are being used to verify compliance with the laser safety standard.
- Coordinating and conducting investigations of laser **incidents** involving serious harm or notifiable events (including serious eye injuries, serious burns, and fires).

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Training and requirements

People can use lasers only after receiving required training (in the case of unrestricted lasers, this may only require the user to read the instruction manual). If people are to use or operate restricted lasers in a workplace, they must be trained in the skills, capabilities and competencies to work safely, as per the *University's Laser Competency Matrix (below)*. Any training received should be recorded on a person's record of learning or similar database.

In the case of an experienced laser operator being recruited to a workplace, they will not require basic laser safety training if they can demonstrate their competency to the laser supervisor. A record of this demonstration of competency must be added to their record of learning.

Types of training within a workplace may include:

- Laser safety awareness training. An e-learning course or short in-house workshop delivered by the LSO in order to communicate University laser safety protocols and to inform participants of basic laser safety information such as mechanisms of laser injury, classes of laser, and basic safety precautions and controls.
- Workplace safety brief: Delivered by the laser supervisor or laser laboratory coordinator, this training can be used to instruct users and operators on correct set-up of a workplace where restricted lasers are being used.
- Specific instruction: Delivered by the laser supervisor or laser laboratory coordinator, this training is delivered to users and operators of specific restricted lasers. It should include working through operating instructions for each different laser that will be used or operated, aligning of the optics and experiments, and enough practical instruction so that the laser can be used safely.

Vision health check

University staff and students who are working in an environment that potentially exposes them to beams of the power of Class 3B and 4 lasers need to have a vision health check. This excludes microscopes and laser systems with the power of the exposed beam within acceptable limits (less than Class 3B). The eye check will be performed by the Auckland Optometry Clinic in Faculty of Medical Science before starting to work and then, biennially. However, in case of presumed pre-existing injury, or upon clinician recommendation, the eye check might have a different frequency. This medical check is free for students and staff (HOD approval is required for staff). Details of the Auckland Optometry Clinic can be found in the "Emergency" chapter.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

University Laser Competency Matrix (v0.1)

Core training element	Performance criteria	User	Operator	Supervisor	Coordinator
1. Demonstrate an	1.1 State how to access University laser	Х	х	х	x
understanding of legislation,	safety information.				
standards and codes of	1.2 Describe the obligations and actions	Х	х	Х	х
practice relevant to laser use	that apply to an individual's specific role or				
within the University.	level of training.				
2. Identify laser hazards.	2.1 Describe the basic anatomy of the eye.	х	x	х	x
	2.2 Describe the properties of laser light.	х	x	х	x
	2.3 Identify classes of laser.	Х	x	Х	x
3. Assess the risks and plan for	3.1 Outline the principles of risk		х	х	x
working with lasers.	management.				
	3.2 Participate in the risk assessment		х	х	x
	process.				
	3.3 Identify sources of information on laser			х	х
	hazards, associated risks, and risk control				
	measures.				
	3.4 Conduct a hazard survey of an area			Х	x
	where restricted lasers are to be used.				

Approved by: Vice-Chancellor

Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer 22/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

	3.5 Conduct a risk assessment of a task or			х	Х
	process where restricted lasers are to be				
	used.				
4. Identify and implement laser	4.1 Identify common laser risk control	x	x	х	x
risk control measures.	measures.				
	4.2 Identify other potential hazards that		x	х	x
	may affect restricted laser use.				
	4.3 Demonstrate the correct selection,		x	х	х
	inspection and fitting of laser safety				
	eyewear.				
	4.4 Demonstrate competence in pre-use		x	х	х
	inspection of restricted lasers and				
	associated equipment.				
	4.5 Demonstrate competence in operating		х	х	x
	restricted lasers and associated equipment.				
	4.6 Demonstrate knowledge of the correct			х	x
	set up and management of laser				
	laboratories.				
	4.7 Demonstrate knowledge of the correct			х	Х
	set up and management of indoor				
	experiments that involve restricted lasers.				

Approved by: Vice-Chancellor

Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer 23/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

	4.8 Demonstrate knowledge of the correct			x	x
	set up and management of outdoor				
	activities that involve restricted lasers.*				
5. Apply first response rescue	5.1 Describe the actions to be carried out in	x	х	x	x
methods.	the event of lasers causing eye injury or				
	suspected eye injury.				
	5.4 Develop emergency first response				x
	plans.				

*4.8. Only required for those people who use restricted lasers outdoors

Recommended Minimum Training:

Unit	Training provider and record	Description	User	Operator	Supervisor	Coordinator
Standard/	keeper					
University						
Workshop						
	Health Safety and Wellbeing staff member	University Laser Safety Awareness*	х	х	x	x
	Supervisor or Coordinator	Laser Lab Safety Brief	х			
	Supervisor or Coordinator	Laser Lab Induction		х	х	х
	Supervisor or Coordinator	Specific competency training in each laser to be used		х	x	x

Approved by: Vice-Chancellor

Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer 24/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

HRHURA	Health Safety and Wellbeing staff member	How to undertake a risk assessment	х	x
HRSORA	Health Safety and Wellbeing staff member	Signing off on risk assessments		x

*Laser safety awareness may take the form of e-learning or a short in-house course.

Records

These records may be paper-based or electronic databases as per local custom. Workplaces where restricted lasers are used should maintain a Laser Operator Register.

Laser Laboratory Authorised Personnel Training and Access Register

Laboratory:

Laboratory Supervisor:

Laser Laboratory Coordinator:

Name	Date Laser Safety Course Completed	Date Vision Health Check	Date on the Job Training Completed	Date Authorised to Operate Laser	Read, Understood and Agree to Comply with SOPs (Signature)	Level of Competency	Laser Laboratory Supervisor (Signature)

Approved by: Vice-Chancellor

Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

User:	Not authorised to use restricted lasers (Class 3B or 4 lasers, or high powered laser pointers) unsupervised.
Operator:	Unsupervised work with restricted lasers permitted, but not permitted to conduct on-the job training (OJT).
Supervisor:	Able to conduct OJT, but does not have authority to deem people competent to use restricted lasers, unless the
	task is delegated by the LLC (not the responsibility).
Laser Laboratory Coordinator:	Able to conduct OJT and deem people competent to use restricted lasers.

26/79

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

PURCHASE IT

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

28/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Safe equipment

If contemplating to purchase a laser make sure that the classification is clear. The lasers need to be classified in accordance with AS/NZS IEC 60825.1.

All lasers must be clearly labelled with a sticker that shows its classification. If a sticker is not present or the reading on the sticker is not clearly readable, the Laser Safety Officer should be contacted to determine whether the laser is safe to use.



Inventory

All restricted lasers must be registered with the Laser Safety Officer.

The LSO maintains a register of all restricted lasers in the University. This ensures that we know where the lasers are being used, and who is able to operate them.

If you are in possession of a restricted laser or you are not sure of the classification contact the Laser Safety Officer.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Appendix 1. - Template to send to the LSO with laser details

	Lasers in use														
Faculty	Department	Building	Room	Туре	Class	Max Output	Pulsed or C/W	Wavelength (nm)		Remarks (microscope, manually incorporated into a microscope, home made)	Commercial or home made	Laser key present?	Contact name	email	
Medical & Health Sciences	Biomedical Imaging Research Unit	502	2802	Diode-pumped Solid State	ЗВ	10mW	сw	567	Melles Griot_85 YCA 010 034	Leica TCS SP2 Microscope	Commercial	Yes			

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Purchase request

All requests to purchase, construct or otherwise obtain restricted lasers must be approved in writing by the Laser Safety Officer.

Some people have placed themselves and others in danger because they have ordered Class 4 lasers directly from the internet without being aware of the dangers involved, so if you need to obtain a restricted laser, you need to contact the Laser Safety Officer in the planning stages.

The LSO can verify that you have the competencies and the facilities to safely use the laser, or they can help you to get the necessary training and resources so that you are safe.

As an example, the template provided in Appendix 1 can be completed and sent via email to the LSO.

Restricted lasers that have not been approved by the Laser Safety Officer must not be used.

Because of the risk of easily inflicted injury, unapproved restricted lasers must not be used. This includes bringing a private high powered laser pointer to an event run by the University. Disciplinary action may be taken against someone who carelessly or recklessly uses a restricted laser.

Lasers must be installed, checked, inspected, cleaned, maintained and adjusted in accordance with manufacturer's instructions.

Reputable laser manufacturers will recommend how to safely install and use a laser, and to keep it operating at peak performance. These recommendations must be followed.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

32/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021



Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

33/79

HSW009.1

34/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

General rules

Supervisors and academic leaders in charge of work areas need to be able to issue instructions and guidance to ensure that everyone is safe. They are allowed to make **reasonable** requests for people to do things, such as:

- Take off specular reflectors such as jewellery or watches
- Put on PPE
- Not touch equipment until given permission to do so
- Turn off cell phones
- Modify their behaviour if it is liable to affect safety

If the requests are not followed, the person in charge can insist that a person is removed from the workplace. A person can also be requested to leave if the person in charge has any other safety concerns.

Note that the requests must be reasonable, so if a person is being asked to do something unsafe (such as being told to remove safety eyewear so they can better observe an experiment), that person can refuse the request without having to fear repercussions.

Restricted laser operations must cease if a supervisor is called away from the workplace for any reason, or if they cannot oversee the work that is being carried out by a person who requires supervision.

People must not use restricted lasers if they are under the influence of, or are impaired by alcohol, recreational drugs or medicines. This includes being influenced by after-affects (such as being hung-over). Medical advice on whether restricted lasers can be operated when taking prescription medicines should be sought on a case-by case basis.

People must not use any laser to dazzle, blind, or distract other people.

This requirement includes those people using unrestricted lasers and low powered laser pointers. There are severe penalties for deliberately dazzling people, especially if they are operating vehicles or aircraft. Offenders can be charged with causing unnecessary danger, and face a term of imprisonment of up to 12 months or a fine of up to \$10,000.

35/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

As far as is reasonably practicable, workplaces in which lasers are used must maintain a satisfactory level of housekeeping.


Bench top clutter must be minimised to reduce the chance of the laser beam being reflected or redirected. Given that many areas where lasers are operated have low ambient light levels, trip hazards must also be eliminated where possible.

Bags should be prohibited in the laboratory. Bags are trip hazards, and may cause people to bend down to get items out. In the process of doing so, a person's head might enter the plane of laser beams.

People using any laser must comply with applicable Safe Work Instructions, and/or other applicable source of information that is needed to perform work safely.

Safe Work Instructions are simple one page reference sheets which will normally be posted near the laser; they list personal protective equipment (PPE), pre and post use checks, and other important information.

Alternately, some complex lasers may have detailed step-by-step operating instructions which must be carefully followed in sequence.



Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Personal Protective Equipment (PPE)

Eyewear is the most common and certainly the most important aspect of personal laser protection, wherever there is some risk of laser exposure. Protective eyewear does not, however, preclude a full safety evaluation and consideration of all alternative means of affording protections - such as total enclosure of the beam, interlocks, beam dumps etc. Laser safety glasses are the last line of defence and not a convenient alternative to avoiding any engineering controls that it may be possible to implement.

Laser safety eyewear incorporates optical filters to reduce the laser radiation to below the level of the maximum permissible exposure, whilst at the same time maintaining a reasonable level of visible light transmission.

The filter damage threshold must be sufficiently high to withstand the maximum possible exposure long enough for appropriate avoiding action to be taken. The absorbance of the filters must be appropriate to the wavelength of the laser(s) used. Clear information must be available to identify the eyewear that is suitable for each laser and whether it is for alignment work or full protection.

Users are responsible for exercising care in the use and storage of eye protection. These items are expensive and easily scratched, degrading performance and comfort. Any defects must be notified to the Laser Laboratory Supervisor requiring the defective item should cease until a safe replacement is made.

Some powerful lasers that are hazardous to skin may require users to wear long clothing. Such clothing should be made of cotton or other flame-resistant material, and must not feature metal or reflective buttons or fasteners.

At a minimum, most restricted lasers will specify that protective eyewear must be worn. The



specifications for wavelength and optical density will be stated in the operating instructions, and users must ensure that the eyewear of the correct specification is

37/79


Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

properly fitted so that there are minimal gaps between the PPE and the face (especially if the user normally wears prescription lenses).

Storage spaces at the point of entry

Storage spaces for specular reflectors such as watches, jewellery, need to be provided at the point of entry so that people do not wear reflective items in the laser controlled area. Phones and other distracting devices can also be left here.



The photograph above shows eyewear storage at entry point. Also note safety folder containing the laser operator register and laser operating instructions.

Warnings!

Alignment goggles may only be used when aligning lasers! Sunglasses *do not* provide laser protection.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Good example of laser eyewear protection





39/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Laser pointers

Laser pointers used at the University within classrooms and lecture theatres must be unrestricted, therefore *clearly labelled* as Class 1 or Class 2, and their power must be less than 1 milliwatt. You do not need to take any special precautions with such laser pointers, other than making sure that you do not deliberately use them to blind, dazzle or distract other people.

Cheap laser pointers without labels must not be used as they may exceed safe power levels.

Laser pointers purchased via the internet

High powered laser pointers can easily cause eye damage, and because irresponsible people have used them to dazzle drivers and pilots, they are now subject to NZ Health Regulations. Severe penalties (including jail) can result if you shine them at vehicles or aircraft, and Australia classifies them as prohibited weapons.

There are very few legitimate reasons to use high powered laser pointers, so any proposal to acquire one must be properly justified. Because they are considered to be restricted lasers by the

University, you must have written permission from the Laser Safety Officer before you buy or otherwise obtain one.

To clarify, a high-powered laser pointer is a device that:

- is commonly known as a laser pointer
- is battery operated
- is designed or intended to be operated while held in the hand
- produces a coherent beam of optical radiation of low divergence
- has a power output of greater than 1 milliwatt



Laser warning label - example

40/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021



5000mW Blue Laser Pointer

Microscopes

Introduction

In terms of laser safety, most microscopes are unrestricted because, although the laser beam is accessible between the objective and the sample stage, the distance between the two is usually very short and crucially, the laser beam is highly divergent and therefore the power at this point is likely to be below the Class 3B level. Generally there are no other points at which the laser beam is accessible.

Summary of hazard

Unrestricted microscope: accessible laser power below the Class 3B minimum Restricted microscope: accessible laser power within the Class 3B or 4 range.

Risk assessment

The Laser Safety Officer will make a risk assessment before you start to work. It may be appropriate to measure the power output at the objective to demonstrate that it is below the Class 3B minimum.

Commercial microscopes

All the commercial microscopes need to undergo an annual service for safety features and accessible laser power by an authorized agent or the manufacturer. The report will be kept in the records by the Laser Laboratory Supervisor.

In-house made microscopes

The Laser Laboratory Supervisor needs to keep records of an annual report containing

- accessible laser power measurements
- check of safety features

The Laser Laboratory Coordinator must ensure that a risk assessment and safe operating procedures are in place. These must take account of all operations - use, cleaning, maintenance and servicing.

During operation

• Do not look directly into the laser beam.

41/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

- Do not remove the specimen while scanning.
- Do not scan without a specimen in position.
- Do not remove any of the optical components including objective lenses without permission from the Laser Laboratory Supervisor
- Do not remove laser keys or attempt to uncouple lasers or safety measures such as interlocks or collars.
- For the **inverted confocal systems**, do not tip the head back when the laser is scanning.
- For the **upright confocal systems**, do not lower the stage to load position or raise the objective to top when laser is scanning.
- When seated, avoid looking continuously into the laser beam (it is scanning at eye level).
- If a sign is on the door indicating that a service engineer is working on the system, **do not enter** before getting permission from the Laser Laboratory Supervisor.
- At end of session, call the next user or shut down if no one is following.
- In the event of exposure, follow the instructions at the "Emergency procedures" section.

Records

- Laser Laboratory Supervisor needs to keep the Risk Assessments and also send them to the LSO
- Laser Laboratory Supervisor needs to keep the reports of the annual services.

Laser survey

For restricted microscopes, the Laser Laboratory Supervisor must carry out an annual Laser laboratory evaluation checklist to check the appropriate control measures are in place (check "University Laser Laboratory Evaluation Checklist").

42/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

BUILD IT/MAINTAIN IT

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

44/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

How to set up a laser lab

The opportunity to set up a new laboratory is always a great occasion, but in doing so, there are items that can be easily overlooked. For this reason, we the Health, Safety and Wellbeing Service have provided a tool that can help University researchers who are undertaking this task. The "How to set up a laser lab" booklet's goal is to identify the do's and don'ts when setting up a laser laboratory. It is available in the Health, Safety and Wellbeing website.

Risk assessment

A risk assessment must be carried out to identify the risks associated with specific restricted lasers.

Laser supervisors and laser laboratory coordinators should have completed risk assessment training so they can identify hazards and controls in the workplace. An example model risk assessment for the operation of a laser laboratory is available as part of the Laser Safety Protocol, and can be amended so that unique circumstances or local conditions can be taken into account. The risk assessment needs to be signed off by a laser laboratory coordinator or the LSO.

Risk assessments should be re-evaluated on the following occasions:

- Annually
- If the workplace undergoes a major change
- If there has been an incident

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Task/Scenario	Associated hazard(s), who may be harmed and how	Existing controls		rent Ratin onsequ Likeliho (R)Rati	g ence x bod	Additional Controls required
			С	L	R	
Untrained or unsafe personnel	 Staff Students Visitors Unsuitable or incompetent operator, unauthorised use. 	 Restricted laser equipment must not be used by personnel who have not passed laser safety awareness training. Operator must be in a fit state to operate the equipment and not impaired by drugs, alcohol or fatigue. Operators are to comply with Safe Work Instructions/Operating Instructions relevant to the equipment. Operators must be supervised where required. 	4	1	4	
Unsafe equipment	 Staff Students Safety critical features 	 All lasers and associated equipment is to be inspected by the operator before the first use each day it is operated. 	4	1	4	

Risk assessment for lasers – screen grab of page 1

Any controls identified in a risk assessment must be implemented so that the residual risk is as low as is reasonably practicable.

You must try to *eliminate* the risks associated with lasers as much as you can, and only adopt a less effective control if it's *unreasonable* to use a higher level control. The ways laser risks can be controlled (in order of most preferred to least preferred method of control) are:

- Elimination totally enclose the laser so that people cannot be exposed to harmful laser energy. This is the most preferred of all the controls and should be used wherever possible.
- Substitution- minimise the risk of harm by using a less powerful laser or less powerful setting to accomplish the task.
- Isolation through engineering minimise the risk of harm by isolating a person from a high risk laser (laser controlled areas, beam guides and curtains).
- Minimisation through administrative controls (procedures and training).
- Minimisation through personal protective equipment (laser safety goggles and appropriate clothing)

For further advice on control measures, contact your supervisor or the LSO.

Any laser-related activities with a residual risk level greater than moderate must not proceed.

This is a fancy way of saying "if it looks unsafe, we need to stop".

46/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

If we are using a restricted laser, we should have first carried out a risk assessment. If any of the controls we identified in the risk assessment are not available, or if circumstances change, we must stop our activity.

Examples are:

- We discover that the interlock to the door of a laser lab has been defective and we can activate the laser with the door open. We must stop the activity until it is fixed.
- Somebody discovers that the provided safety eyewear does not fit over their prescription lenses. The activity must not start until they are properly protected.
- A laser light show has been set up so that the lasers project three metres above the audience. We notice that a woman is sitting on her boyfriend's shoulders and is holding a selfie stick. The laser show must be stopped and she must be at ground level before we can continue.

Example risk assessment

Example risk assessments can be found on the Health, Safety and Wellbeing website.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Health, Safety and Wellbein	ng (HSW) Risk Assessment		
Document No:		Assessment Date:	
Faculty/ Service Division:		School/Department :	
HSW Risk:	Uncontrolled Risk is High, Controll	ed is Moderate	
Form completed by:		Responsible Line Manager:	
Signed		Signed	
Dated		Dated	
Other Risk Assessments which might also be required:	Manual Handling	Working at Night	Workstation
Description of activity and/or location:	Use of Restricted Lasers in Laborat	tories	
	Potential Hazards – Beam related haza Potential Harm – Blindness, skin burns		

Identify Hazards and Control the Risks:

1. An activity may be divided into tasks. For each task identify the hazards and associated risks. Also list the possible scenarios which could sooner or later cause harm.

2. Determine controls necessary based on University standards, legislation, codes of practice, AS / NZ standards, manufacturer's instructions etc.

3. List existing risk controls (take credit for what you do)

4. Rate the risk once all controls are in place using the matrix in

5. List any additional controls that need to be implemented and take action

6. Communicate the findings

The boxes will resize to suit your situation/the amount of text you need to use - press tab after last cell to create new rows

Task/Scenario	Associated hazard(s), who may be harmed and how	Existing controls	Risk Ra (C)Consequ (L)Likelih		Current Risk Rating (C)Consequence x (L)Likelihood = (R)Rating C L R		ting ence x pod ng	Additional Controls required
Untrained or unsafe personnel	 Staff Students Visitors Unsuitable or incompetent operator, unauthorised use. 	 Restricted laser equipment must not be used by personnel who have not passed laser safety awareness training. Operator must be in a fit state to operate the equipment and not impaired by drugs, alcohol or fatigue. Operators are to comply with Safe Work Instructions/Operating Instructions relevant to the equipment. Operators must be supervised where required. 						
Unsafe equipment	StaffStudents	• All lasers and associated equipment is to be inspected by						

	Safety critical features may be missing or damaged.	 the operator before the first use each day it is operated. Non-conformities to be reported to supervisor and rectified before being used. 		
Inappropriate facilities	 Staff Students Inappropriate facilities may not contain laser beam hazards. 	 All use of restricted lasers must take place in a designated laser laboratory that has been evaluated by the Laser laboratory Coordinator. Laboratories must be suitable for the lasers being used. 		 The Laser Supervisor is legally responsible for the work being undertaken. Supervisors should ensure that laboratory infrastructure is safe for the lasers being used (eg. laser controlled areas, non-flammable non-reflective surfaces, curtains, safety devices, etc).
Personal Injury	 Staff Students Laser energy may penetrate the eye or burn skin. 	 All personnel observing the laser must wear protective eyewear. All personnel must remove personal specular reflectors upon entry to the laboratory. Where required, non-flammable long clothing must be worn. Laser beams must not be placed at eye level (sitting or standing). Beam guides and dumps should be used where possible. Laser arming keys must be kept under positive control at all times by the person in charge of the laser. 		 The laser operating instructions will specify the eyewear requirements for optical density and wavelength. Eyewear must fit the user, especially if the user wears prescription lenses. All clothing must be free of specular reflectors such as buttons and nametags. Laser arming keys must only be inserted just prior to the laser being activated. Lone work with restricted lasers is not permitted.

Fumes and vapours.	 Staff Students Toxic fumes can cause long term harm. 	 Ensure that personnel are aware of the hazards associated with the experiment. Do not create fumes unless it is part of the experiment. Arrange the work so that fumes do not rise into the worker's face. Use local exhaust ventilation. Use personal protective equipment (PPE). 		• Any respiratory protection must be appropriate for the fumes being created.
Fire	 Staff Students Laser beams may ignite combustible materials. 	 All combustible materials to be removed or made safe. No flammable liquids, vapours, gases or dusts to be present. Extinguishers/hoses to be provided outside the lab. Operator must know how to use fire equipment. Operator must know location of telephone/fire alarm. Site to be inspected after completion of work for smouldering surfaces or objects. 		
Electrical	 Staff Students High Voltage equipment may deliver electric shocks. 	 All lasers and associated electrical equipment to be checked in accordance with the University's electrical appliance testing protocol. Non-conformities to be reported to supervisor and rectified before being used. 		

Action Plan						
Management agreed	Resources	Action By:		Action Complete: Responsible Line Manager		
additional control measures to be implemented	Required	Responsible Person	Target Date	Target Date Completion Date		Date

Review	
Review Details	Comments
Scheduled Review Date	
Are all control measures in place?	
Are controls eliminating or minimising the risk?	
Are there any new problems with the risk?	
Are the supervisory arrangements adequate?	
Are the levels of skills, capabilities and training	
adequate?	
Review By: (name)	
Review Date:	

Communication				
	Method	Yes	Date	Comments
Reference of formal communication to staff	Copy of risk assessment issued to staff			
	Controls covered in team procedure issued to staff			
	Staff handbook issued to staff			
	Other			
How they were consulted on the risk	Health, Safety and Wellbeing Committees			
Additional Methods of	Induction			
Communication	Toolbox Talk			
	Team Meeting			
	Email circulation			
	Other			

115 0	N Risk Assess	sment Matrix				
	4 Probably	ery likely expect the event to most circumstances	Moderate (4)	High (8)	Extreme (12)	Extreme (16)
Likelihood Level		Likely ely to occur at least er the coming year	Moderate (3)	High (6)	High (9)	Extreme (12)
(elih	2	Possible	Low	Moderate	High	High
5		occur at some time	(2)	(4)	(6)	(8)
	1 Occurrei	Unlikely nce is conceivable, expected to occur	Low (1)	Low (2)	Moderate (3)	Moderate (4)
			Minor	Moderate	Major	Severe
			1	2	3	4
				Conseque	nce level	
Consequence description	Harm to Pe Potential for	ople injury or death	None or trivial / negligible injury (no or slight injury which requires localised first aid)	Conseque Minor injury (illness or injury is not serious, medical treatment required)	Serious injury (serious injury or illness, hospitalisation required)	Fatality, major injury (death, permanent disablement, or significant long-term illness)

_	None or issue raised by	Internal scrutiny to	Medium-term	Persistent stakeholder
Reputation and Legal		,		
Potential for publicity with a	staff or students and	prevent escalation and	stakeholder concern,	concerns, international
negative impact on	resolved promptly by	short-term stakeholder	national media	media scrutiny and long
reputation / potential for	management	concern	scrutiny and 'brand'	term 'brand' impact
legal prosecution			impact	
	None or legal dispute –	Minor non-compliance,		Significant non-
	found not guilty – fines	limited notification to	Medium non-	compliance, extensive
	up to \$3 million (Body	regulators / affected	compliance, moderate	notification to regulators
	Corporate), \$600,000	stakeholders	notification to	/ affected stakeholders,
	(Officer)		regulators / affected	potential for legal
			stakeholder, potential	proceedings /
			for legal proceedings /	imprisonment / fines
			fines	
			-	
Operations	None or business	Business interruption		
Extent of ability to maintain	interruption < 4 hours	between 4 hours to 5	> 5 days	many weeks
core business		days		
	None or effectiveness	Operational disruption	Medium operational	Breakdown of key
	and efficiency of a	manageable by	impact resulting in	activities and significant
	service, programme or	workarounds	delay of key	long-term impact
	project impacted in the		deliverables	
	short term			
		Moderate damage to		Massive damage to
	None or slight damage to	property or equipment	Major damage to	property or equipment
	property or equipment		property or equipment	
			Field of addition	

Environment	None or minimal impact	Minor short-term or	Serious, medium-term	Very serious, long-term
Extent of negative impacts on the environment		intermittent impact, able to be contained with specialist		or permanent damage
		assistance		
	None or clean up		Clean up expenses up	Clean up expenses >
	expenses up to \$25,000	Clean up expenses up	between \$1m - \$5m	\$5m
		between \$25,000 to		
		\$1m		

Consider the Likelihood

How often is the task done? Has an accident happened before (here or at another workplace)? How long are people exposed? How effective are the control measures? Does the environment affect it (e.g. light, temperature, space)? What are people's behaviours (e.g. stress, panic, deadlines)? What people are exposed (e.g. disabled, young students, etc.)?

Consider the Consequences

What type of harm could occur (minor, serious, death)? Is there anything that will influence the severity (e.g. proximity to hazard, person involved in task, etc.)? How many people are exposed to the hazard? Could one failure lead to other failures? Could a small event escalate?

Calculate the Risk

The final score for each risk is calculated by multiplying the likelihood and consequences response scores. This will give a risk score of between 1 and 16.

All risks rates as "High" or "Extreme" require detailed analysis of mitigating practices / controls to determine the residual risk rating. Action must be taken.

"Low" and "Moderate" risks may be excluded from further analysis (other than when the consequence may be severe). However the rationale for excluding these risks should be documented to demonstrate the completeness of analysis undertaken. **Some action may be required.**

Other than in the most unlikely circumstance, risks that can cause major or severe harm to people have been determined as "high" or "extreme". Management review is considered appropriate for risks of these nature due to the potential magnitude of the impact, even though the likelihood may be assessed as relatively low.

HSW Risk Assessment Form

Risk Priority - Legend

Extreme	Intolerable risk. Immediate action(s) is to be taken by Faculty/Service HSW risk owners - including DVCs, Deans of Faculties,
(12-16)	Directors of Services, Academic Heads/PIs, Services Managers. Work should not be started or continued until the risk has been
	reduced to as low as reasonably practicable using the hierarchy of risk controls. The Associate Director Health, Safety and
	Wellbeing, and Manager Risk and Performance must be advised of the risk for their review. The risk should be included in the
	UoA wide risk register.
High	Should not be tolerated. Urgent action is to be taken by the immediate manager. Work should not be started or continued
(6-9)	until the risk has been reduced to as low as reasonably practicable using the hierarchy of risk controls. The HSW Manager
	working with the Faculty/Service, and Manager Risk and Performance must be advised of the risk for their review. To be included
	in the UoA wide risk register.
Moderate	Management to monitor risks in case changing circumstances increase the level of risk. Some action may be required, e.g.
(3-4)	improving controls.
Low	Requires no further attention above routine practices and procedures, apart from monitoring.
(1-2)	

Note: This proposed Health and Safety Risk Assessment Matrix aligns with WorkSafe NZ guidance, UoA Resilience Management Plan, UoA Risk Determination Matrix, UoA TVRA and UoA Incident Levels

Repairs and maintenance

Records

Where lasers and associated equipment require recorded scheduled servicing, maintenance activities must be recorded. If a laser requires a repair, we need to know who repaired it and what work was done to the laser to fix it. All servicing must be recorded by the laser laboratory coordinator or the laser supervisor (if delegated by the LLC).

Class 3B and 4 laser systems must undergo a check, before being put into use for the first time and annually thereafter. A template can be found in the "University Laser Laboratory Evaluation Checklist" chapter.

Temporary laser controlled areas

A temporary laser controlled area should be established whenever conditions allowing human access to hazardous levels of laser radiation are created temporarily, (e.g. during servicing), and where persons who are unauthorized, unaware of the presence of the laser hazard and/or are not appropriately trained or supervised in the necessary safety procedures could be present.

The guidance for temporary controlled areas are the same as for the laser controlled areas (refer to "How to set up a laser lab"). Although the normal requirement for engineering control of access may be difficult to achieve, administrative controls can have increased effectiveness when restriction of access is only temporary. If safe access is not controlled by engineering means, then appropriate warning and prohibited entry signs should be posted at the point of entry.

Controls during servicing

- Comfortable laser safety goggles with correct wavelengths for use by Laser Laboratory Supervisor or Coordinator to perform repairs or maintenance.
- Reducing the level of emission to the maximum necessary.
- Limiting the range of movement of beam steering components.
- Placing large area beam stops behind the target to stop the laser beam in case of missed target.
- Using non-reflective coating on tools and requiring the removal or covering of jewellery or watches.

59/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

60/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021



Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

61/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

62/79

Physical security and access

All laboratories and storage areas containing restricted lasers in the University of Auckland must have access control. At a minimum this will mean locked doors which are opened using conventional keys. However, swipe card access is preferable. Where segregated swipe-controlled access is not available, areas should be controlled via key or combination locks.

Information security

Why do we need information security?

The University has a duty of care to protect the privacy of University of Auckland staff, students and stakeholders (such as contractors), as well as to ensure confidential information is only viewed by those authorised to do so. Information security is also designed to ensure that the data is kept intact and accurate for periods of up to 10 years.

What information needs to be kept securely?

The requirements for information security outlined in this document apply to databases or protocols which contain:

- • Personal identities of staff and students at the University.
- • Information about access control and contingency plans for the facilities.
- • Confidential memorandums and emails.

This includes (but is not limited to):

- • Training databases.
- • Restricted lasers databases.
- • Audits databases and corrective actions.
- • Containment manuals (including archived versions).

Where is the information held?

All information related to University of Auckland laser laboratories is to be held on information technology equipment that is:

• Exclusively owned and controlled by the University of Auckland.

63/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

- Only accessible using University of Auckland-controlled networks.
- Only accessible using University of Auckland passwords.
- Regularly backed up onto University-owned and controlled drives.

Using shared databases

Where more than one person requires access to this information, it is to be held on separate password-controlled shared drives.

Shared databases are to be configured in a robust manner to prevent corruption of files. Wherever possible SQRL databases and web-based access is to be used.

If external virtual drives (e.g. DropBox or data storage in the Cloud) are used, the security of data storage is to meet the security requirements of the University of Auckland Information Security services.

Review of information security

Information security is to be reviewed annually by the University of Auckland Information Security services.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

MOVE IT/DISPOSE OF IT

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

65/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

66/79

Transfer

If you want to relocate a functional restricted laser, or it is no longer needed and you would like to donate it, you need to do so through a written communication to the LSO (an email suffices).

You can transfer or sell the laser within the University or to a responsible organisation (such as another university). Note that the laser system must be serviceable and fit for use. It should be sent with all instruction manuals and other relevant user information. In addition, you must check that the receiver has a viable laser safety programme and a laser safety officer in place before you transfer the laser. *Restricted lasers are not to go to public auction*.

The LSO will approve the transfer and will update the restricted lasers register.

Disposal

Defective lasers and associated equipment are to be taken out of service as soon as defects are identified. This is especially important if there is any damage or defect that affects the ability to align the laser, or if there is damage to any electrical insulation or power cords. Lasers being decommissioned must be disposed of responsibly. All lasers are considered to be electronic waste, and unrestricted lasers may be taken to an e-waste recycling centre once any batteries have been removed.

When restricted lasers are being disposed of the Laser Safety Officer must be notified. This will ensure that the laser register is up to date, and that the selected disposal method is appropriate.

Options for the disposal of restricted lasers are:

- Send it back to the manufacturer for recycling or disposal. Not all manufacturers will do this, and there may be a cost involved.
- Dispose of the laser. If this option is selected, the laser must be sent to the Laser Safety Officer for destruction.

Any costs involved will need to be borne by the unit disposing of the laser.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

68/79

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

MEASURE IT

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

69/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

70/79

Internal verification

Regular internal **verification** is an integral part of any risk management system. The purpose is two-fold:

1. To ensure that systems and procedures outlined in containment guidance are fit for purpose and achieve the desired results; and if not, appropriate improvements are identified

2. To identify risk and ensure appropriate risk mitigation procedures are put in place

Verification consists of:

- Annual Laser laboratory inspection by Laser Laboratory Coordinator.
- Annual verification of the restricted lasers registered by the LSO.
- Annual verification of records of trained users by LSO.
- Annual verification of the annual safety service for microscopes by LSO.
- Unannounced laser laboratory checks.
- Identification of underlying issues in order to ensure the reasons are addressed.
- Identification of exemplars (practices or facilities).

The results of verification exercises are as follows:

- Where **non-conformances** are identified they are documented and corrected.
- The root cause of non-conformance is also identified (where appropriate), documented and corrected.
- Exemplars identified are noted in the facility's verification register and used as models for improvement (i.e. good practices are shared across the University community).

Annual laser laboratory inspection

Annual Laser laboratory inspections are to be undertaken by Laser Laboratory Coordinator or Supervisor. These are designed to ensure that the physical structure of laboratories is maintained to a high standard, that the safety features are in place and that individuals are routinely observing good laboratory practice. A checklist is provided in "University laser laboratory evaluation checklist".

Annual verification of the restricted lasers register

Restricted laser, commercial or home-made systems, must be registered with the LSO. The LSO will arrange a meeting every year with the Laser Laboratory Coordinators to verify the accuracy of the Restricted Lasers Register. During the meeting, LSO will also check the trained users' records and the annual safety service for microscopes.

Documentation in each laser laboratory

Each laboratory is to have the following registers as hard copy:

- University laser laboratory evaluation checklist.
- Laser Laboratory Authorised Personnel Training and Access Register.
- List of the restricted lasers.
- Annual safety service (only for microscopes).
- Risk assessment (if any).

University laser laboratory evaluation checklist

To be completed for Class 3B and 4 Laser Systems, before being put into use for the first time and annually thereafter.

Location:

Dated:

Element	Y/N/NA	Remarks
Laser Identification*		
Description		
Class		
Wavelength		
Maximum Output		
NOHD (if known)		
All restricted lasers registered with		
the Laser Safety Officer		
Administration	· · · ·	
Laser operator register displayed		

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing

Content Manager: Laser Safety Officer

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Laser supervisor(s) identified				
UoA laser safety standard displayed				
Laser activation flowchart displayed				
Operating instructions				
available/current				
Spectators briefed and supervised				
Entrance to Laser Laboratory				
Appropriate warning signs				
Entry way controls				
• Interlocked doors (where poss	ible)			
• door/curtain at entry way				
• Visible/audible signals (e.g. il sign or warning lights)	lluminated			
An area for protective eyewear and				
specular reflectors provided near				
entrance				
Appropriate eye protection clearly				
identified				
Skin protection when needed (i.e.				
UV lasers)				
Windows/doorways				
covered/restricted				
Panic button/e-stop available				
Laser controlled area within a				
laboratory				
Laser hazard area clearly				
identified				
Laser beam path enclosed				
Laser and Optical Components				
Laser warning signs and labels				
legible				
Protective housing and shields intact				

73/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Keys under positive control				
Beam stop/attenuator provided				
Reflective materials out of paths way				
Laser secured to table or base				
Optical components secured to table				
or base				
Beam routed above or below eye				
level				
Non-Beam Hazards				
Electric shock/high voltage				
minimised				
Fire potential minimised				
Vapour/particle generation				
minimised				
Hazardous or toxic chemicals				
minimised				
Metal laser tables earthed				
Check carried out by	Signature		Date	

*Only restricted lasers need to be identified and registered.

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System



Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

75/79

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

76/79

Emergency procedures

The University has published procedures for common emergencies such as fire, earthquake and power cuts, however a laser lab must have extra procedures in place in case someone is exposed to harmful laser radiation. A simple laser emergency plan has been developed as an example:

- Stop the beam with the emergency button
- Tell your workmate you have been exposed to the laser and have hurt your eyes or skin
- Get them to help you or call Security at ext 85000. Ask to be escorted to the Auckland Optometry Clinic in the Faculty of Medicine as soon as possible
 - Auckland Optometry Clinic operating hours: Monday to Friday 8am-5pm (Phone 09-923-9909)
 - Out of hours: go to the Emergency Department at the Auckland Hospital
- Do not disturb the equipment setup
- Cover your injured eye/s before going into daylight

After the accident complete an Incident Form and inform your supervisor and LSO.

Master stop switches (laser kill switches) should also be installed at strategic positions throughout the laser lab.

Incident reports

Workplaces must record and report incidents (including accidents or near-misses) to ensure that any learning from such an incident is captured and not repeated. Incidents that cause serious harm or are notifiable incidents must be reported to the HSW Service without delay, and the accident scene must be preserved as much as possible for an investigation.

In the case of laser injuries, the following are notifiable to WorkSafe NZ:

A serious eye injury that requires	Includes:
immediate treatment (other than first	an injury that results in, or is likely to
aid)	result in, the loss of an eye or vision
	(total or partial)

77/79

Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

	Does not include: exposure to a substance or object that only causes discomfort to the eye
A serious burn that requires immediate treatment (other than first aid)	Includes: A burn that needs intensive or critical care such as a compression garment or skin graft.
	Does not include: A burn treatable by washing the wound and applying a dressing.

A notifiable incident is where someone's health or safety is seriously endangered or threatened. People may be put at serious risk even if they were some distance from the incident. Accidentally lasing a member of the public when using a high powered laser pointer outdoors could be considered a notifiable incident.


If there is any doubt about whether the incident needs to be notified, contact the Health, Safety and Wellbeing Service on 09-9234896 as soon as possible. Ask for the Laser Safety Officer.

Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System

Appendix 1. How to find the Eye Clinic

The Auckland Optometry Clinic is located at 143-151 Park Rd, Grafton, Auckland 1023, Building 505, Level 3.



Approved by: Vice-Chancellor Document Owner: Associate Director, Health, Safety and Wellbeing Content Manager: Laser Safety Officer Version: 1 Issue Date: 19/02/2018 Review Date: 19/02/2021

Once printed this document is uncontrolled. Health Safety and Wellbeing Management System